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Droga Test Cocaina (COC)

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Result in 5 minutes

INTENDED USE AND SUMMARY

Cocaine (COC) Rapid Test is a lateral flow chromatographic immunoassay for the detection of the cocaine metabolite, benzoylecgonine in human urine at a cut-off concentration of 300 ng/mL. This test will detect other related compounds, refer to the analytical specificity table in this insert.

This test provides only a preliminary result of the analytical test. It is necessary to use a more specific alternative chemical method to obtain a confirmed analytical result. Gas chromatography and mass spectrometry (GC/MS) is the method of confirmation to be preferred. Any findings for the determination of the presence of drug abuse should be related to clinical considerations and professional judgment, especially when the preliminary result is positive.

Cocaine is a powerful central nervous system stimulant (CNS) and a local anesthetic. Initially it causes extreme energy and restlessness while gradually results in tremors, hypersensitivity and spasms. In large quantities, cocaine causes fever, numbness, difficulty breathing and loss of consciousness. Cocaine is often self-administered for nasal inhalation, intravenous injection, and free smoking. It is excreted in the urine in a short time mainly as benzoylecgonine. Benzoylecgonine, an important metabolite of cocaine, has a longer biological half-life (5-8 hours) than cocaine (0.5-1.5 hours) and can generally be detected for 24-48 hours after exposure to cocaine.

The rapid test for cocaine (COC) is a urine screening test that can be performed without the use of special tools. The test uses a monoclonal antibody to selectively detect high levels of cocaine metabolites in the urine if it exceeds the cut-off level 300 ng/ml.

This is the screening limit suggested for positive samples established by the Substance Abuse and Mental Health Service Administration (SAMHSA, USA).

PRINCIPLE

The Rapid Test for Cocaine Detection (COC) is a rapid immunochromatographic test based on the principle of competitive binding. Drugs that may be present in the urine sample compete with the drug conjugate for antibody binding sites.

During the test, a sample of urine migrates upwards by capillary action. Benzoylecgonine, if present in the urine sample below 300ng/ml, will not saturate the sites: saturation of the binding sites of antibody coated particles in the test strip. Antibody-coated particles will then be captured by the immobilized benzoylecgonine conjugate and a visible colored line will appear in the test line (T) area. The colored line will not form in zone T if the benzoylecgonine level exceeds 300ng/ml because it will saturate all antibody binding sites against benzoylecgonine.

A drug-positive urine sample will not generate a colored line in zone T due to drug competition, while a drug-negative urine sample or a sample containing a drug concentration below cut-off will cause the formation of the colored band in the specific test area. As a control of the procedure, a colored band will appear in the reactive zone, indicating that a correct amount of sample has been used and that the membrane migration has occurred.

COMPOSITION

Test Disposable cocaine, instructions for use, jar for collecting samples.

STORAGE AND HANDLING

- Store the test in a cool, dry place at a temperature between 2 and 30° C. Keep away from light. Exposure to temperature and/or humidity outside the specified conditions may cause inaccurate results.
- Do not freeze. Use the test kit at temperatures between 15 and 30° C.
- Use the test kit between 10 and 90% humidity.
- Do not use the test kit after the expiration date (printed on the bag and package).

N.B. All expiration dates are printed in Day-Month-Year format.

WARNINGS PRECAUTIONS AND LIMITATIONS

- Product not intended for medical or diagnostic use.
- Do not use after the expiration date.
- Keep the panel of test strips in the casing until ready for use.
- All samples must be considered potentially dangerous and therefore must be handled with the precautions of use relating to potentially infectious products.

- This kit contains products of animal origin. The certified knowledge of the origin and/or health status of the animals does not fully guarantee the absence of transmissible pathogens. It is therefore recommended to treat these products as potentially infected and to handle them according to the usual safety precautions (do not ingest or inhale).
- Cross-contamination of samples shall be avoided by using a new sample collection container for each specimen obtained.
- Carefully read the entire procedure before performing any tests.
- Do not eat, drink or smoke in the area where samples and kits are handled.
- Observe established precautions against microbiological hazards throughout the procedure and follow the standard procedures for proper sample disposal.
- Wear protective clothing such as lab coats, disposable gloves and eye protection during sample analysis.
- Humidity and temperature can adversely affect the results.
- The materials used for testing must be manufactured in accordance with local, state and/or federal regulations.
- It is possible that technical or procedural errors, as well as other interfering substances in the urine can cause incorrect results.
- Adulterants, such as bleach and/or alum, may produce incorrect results in urine samples regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine sample.
- A positive result indicates the presence of the drug or its metabolites but does not indicate the level of intoxication, route of administration or concentration in the urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when the drug is present but below the cut-off level of the test.
- The test is not able to distinguish between the drug and a medicine containing the same substance, false positives can be caused by some drugs.

SPECIMEN COLLECTION AND PREPARATION

1) Analysis of urine

The urine sample must be collected in the appropriate jar inside the clean and dry package. Urine can be collected at any time of the day. Urine samples with visible precipitates should be centrifuged, filtered or allowed to settle to obtain a clear sample for testing.


2) Storage of samples

Urine samples can be stored at 2-8 C up to 48 hours before analysis. For prolonged storage, samples may be frozen at a temperature below -20°C. Frozen samples must be thawed and mixed before testing.


TEST PROCEDURE

Allow the panel and urine sample to reach room temperature (15-30°C) before testing.

- Remove the test device from the sealed bag and use it as soon as possible.
- Place the test device on a flat surface and drive. Hold the dropper vertically and transfer 3 whole drops of urine (about 120µL) into the sample well (S) of the test device, then start the timer. Avoid trapping air bubbles inside the cockpit. See illustration below.
- Wait for the red lines to appear. The result should be read after 5 minutes. Do not interpret the result after 10 minutes.

INTERPRETATION OF TEST RESULTS

POSITIVE: A colored line in the control line area (C) but no line in the test area (T) for a specific drug indicates a positive result. This indicates that the drug concentration in the sample exceeds the designated cut-off for that specific drug.

NEGATIVE: Two distinct colored lines appear. A colored line in the control line area (C) and a colored line in the test the line region (T) for a specific stave indicates a negative result. This indicates that the drug concentration in the sample is below the threshold level designated for that specific drug.

N.B. The color gradient in the reactive zone (T) may vary, but must be considered negative whenever there is also a weak colored line.

INVALID: The control line (C) is not colored. The most likely reasons are insufficient sample volume or incorrect

procedural techniques. Review the procedure and repeat the test using a new panel. If the problem persists, immediately stop using the test and contact your local distributor.

QUALITY CONTROL

The Test includes an internal control system: a colored line appears in the control zone (C). This coloured band indicates that sufficient flow has occurred and the functional integrity of the test device has been maintained. The control standards are not supplied with the kit; however, it is recommended that negative and positive controls be tested by a laboratory to confirm the test procedure and verify its performance.

Produced and distributed by:
Brands2China S.r.l.
brand owner:
Protezionipiù
Via Roma, 33
36014 - Santorso - Italia
www.protezionipiu.com