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Multi-droga Test Urina 5

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Multi-drug Test Urine 5

To perform with urine Result in 5 minutes

A rapid, one step screen test for simultaneous, qualitative detection of multiple drugs and metabolites in human urine: amphetamine, methamphetamine, cocaine, marijuana and opiates.

WARNINGS PRECAUTIONS AND LIMITATIONS

- Not for medical or diagnostic use.
- Do not use after expiration date.
- Test device should remain sealed until ready for use.
- All specimens must be considered potentially infective products.
- Used testing materials should be discarded in accordance with local regulations.
- The drug test provides only a preliminary analytical result. A more specific chemical method must be used.
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleach and/or alum, in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- This test provides only preliminary analytical data. More specific chemical methods should be used as confirmation of the analytical result. Gas chromatography and mass spectrometry (GC/MS) are considered preferred confirmatory methods.
- Clinical considerations and professional judgment should be applied to any drug of abuse result, particularly when the preliminary result is positive.
- The test does not distinguish between drugs of abuse and certain medications.
- A positive result might be obtained from certain foods or food supplements.

INTENDED USE AND SUMMARY

The Rapid Multi-Drug Disposable Test is a lateral flow chromatographic immunoassay for the qualitative detection of the following drugs without the aid of instruments.

Test	Calibrator	Cut-off (ng/ml)
Cocaine(COC)	Benzoyllecgonine	300
Amphetamine(AMP)	d-Amphetamine	1000
Methamphetamine(MET)	d-methamphetamine	1000
Marijuana(THC)	11-nor- Δ^9 -THC-9 COOH	50
Opiates(OPI)	Morphine	2000

PRINCIPLE

The Multi-Drugs Rapid Test is an immunoassay based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody.

During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody coated on the particles. The antibody coated particles will then be captured by the immobilized drug conjugate and a visible colored line will show up in the test line region of the specific drug strip. The colored line will not form in the test line region if the drug level is above its cut-off concentration because it will saturate all the binding sites of the antibody coated on the particles.

A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

COMPOSITION

Urine 5 Multi-Drug Test, Instructions for Use, Sample Collection Jar.

STORAGE AND HANDLING

- Store the test in a cool, dry place at a temperature between 2 and 30° C. Keep away from light. Exposure to temperature and/or humidity outside the specified conditions may cause inaccurate results.

- Do not freeze. Use the test kit at temperatures between 15 and 30° C.
- Use the test kit between 10 and 90% humidity.
- Do not use the test kit beyond the expiration date (printed on the pouch and box).

NOTE: All expiration dates are printed in Day-Month-Year format.

SPECIMEN COLLECTION AND PREPARATION

1) Urine Assay

The urine specimen must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear supernatant for testing.


2) Specimen Storage

Urine specimens may be stored at 2-8°C for up to 48 hours prior to assay. For prolonged storage, specimens may be frozen and stored below -20°C. Frozen specimens should be thawed and mixed before testing.


TEST PROCEDURE

Allow the test, urine specimen, and/or controls to equilibrate to room temperature (15-30°C) prior to testing.

1) Remove the test panel from the sealed bag and use it as soon as possible.

2) Take off the cap outside of the test end. With arrows pointing toward the urine specimen, immerse the test panel vertically into the urine specimen for at least 10-15 seconds. Immerse the test panel to at least the level of the wavy lines on the strip(s), do not pass the arrows on the test panel when immersing the panel. See the illustration below.

3) Put the cap back on the cassette and place the test panel on a non-absorbent flat surface, start the timer and wait for the colored line(s) to appear.

4) Read the results at 5 minutes. DO NOT INTERPRET RESULT AFTER 10 MINUTES.


INTERPRETATION OF TEST RESULTS

POSITIVE: A colored line in the control line region (C) but no line in the test line region (T) for a specific drug indicates a positive result. This indicates that the drug concentration in the specimen exceeds the designated cut-off for that specific drug.

NEGATIVE: Two distinct colored lines appear. A colored line in the control line region (C) and a colored line in the test line region (T) for a specific drug indicate a negative result. This indicates that the drug concentration in the specimen is below the designated cut-off level for that specific drug.

Note: The shade of color in the test region (T) may vary, but it should be considered negative whenever there is even a faint colored line.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test panel. If the problem persists, discontinue using the lot immediately and contact your local distributor.

QUALITY CONTROL

Internal procedural controls are included in the test. A colored line appearing in the control region (C) is the internal procedural control. This procedural control line indicates that sufficient flow has occurred, and the functional integrity of the test device has been maintained. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance.

Produced and distributed by:
Brands2China S.r.l.
brand owner:
Protezionipiù
Via Roma, 33
36014 - Santorso - Italia
www.protezionipiu.com